## Appendix How to Find your Ethernet Address


Several GE Healthcare software products present the Ethernet address during installation, start up or by using a separate tool. Here is a general description on how to find the ethernet address without using any of these tools.

**Important.** Use the Ethernet address for the computer where the license file will be placed. The placement may depend on the license type:


- 1. Concurrent license (a.k.a. floating), use Ethernet address for the computer where the license server is installed.
- 2. Node locked license, use Ethernet address for the computer where the software application is installed.

## **Get Ethernet address**

1. Start->Run



2. Type cmd and click OK



3. In the Command Window type "ipconfig /all" and press enter.

```
C:WINDOWS\system32\cmd.exe

Microsoft Windows XP [Uersion 5.1.2690]
(C) Conveight 1985-2001 Microsoft Corp.

U:\ipconfig /all_
```

4. The response should contain several parameters for the network. Look for the line **Ethernet** adapter Local Area Connection. Then look for the **Physical Address.** This number is the Ethernet address. In the example below it is marked in white. It must be used as a straight number without "-". In the example below it is **00aa11bb22cc.**